



Runtime Verification for Interconnected Medical Devices

Martin Leucker Malte Schmitz Danilo à Tellinghusen

{leucker, schmitz, tellinghusen}@isp.uni-luebeck.de

Institute for Software Engineering and Programming Languages,
University of Lübeck, Germany

7th International Symposium on Leveraging Applications of
Formal Methods, Verification and Validation

Outline

1. Interconnection of Medical Devices
2. Runtime Verification for Medical Devices
3. **Device Modeler**
Risk Analysis through Contract Enforcement
4. **Swiss Army Knife**
Interconnection Debugging Tool

Operation Room 1956



Bundesarchiv, Bild 183-35596-0002
Foto: Höhne, Erich; Polli, Erich | 17. Januar 1956

Operation Room 2016




Safe Interconnection of Medical Devices

Upcoming IEEE 11073-SDC standard

- ▶ Interconnect medical devices **in the operation room.**
- ▶ Dynamic interconnection of devices **from different manufacturers.**
- ▶ Devices announce themselves in network **with an interface description.**




Open Surgical Communication Protocol (OSCP)



Source: Martin Kasparick, Stefan Schlichting, Frank Golatowski, Dirk Timmermann:
New IEEE 11073 standards for interoperable, networked point-of-care Medical Devices. EMBC 2015

Semantic Interoperable Domain Information & Service Model



Source: Martin Kasparick, Stefan Schlichting, Frank Golasowski, Dirk Timmermann:
 New IEEE 11073 standards for interoperable, networked point-of-care Medical Devices. EMBC 2015

Risk Analysis

- ▶ European Medical Device Directive demands **execution of risk management.**
- ▶ How to do risk management for **dynamic interconnection?**



Risk Analysis: Contract Enforcement

- ▶ European Medical Device Directive demands **execution of risk management**.
- ▶ How to do risk management for **dynamic interconnection**?


Embedded **monitors check at runtime** that

- ▶ the **device itself** and
- ▶ **other devices** controlling it

satisfy the **interface description**.



Our Tools



OSCP Swiss Army Knife

UNIVERSITÄT ZU LÜBECK
 INSTITUTE FOR SOFTWARE ENGINEERING
 AND PROGRAMMING LANGUAGES

OSCP Swiss Army Knife

Manage Monitors

Search Devices


Show Monitor Events

Element	Term Code	Info	Attribute	Value
Network			Element	Numeric Metric
MDS	MDC_DEV_ANALY_SAT...	urn:c...	Term Code	MDC_PULS_OXIM_SAT...
Raw MDDescription			Handle	spo2
Activate Operation D...	MDCX_START_IDENTIF...		Observation Time	2016-05-02T23:00:55.301
Activate Operation D...	MDCX_STOP_IDENTIFIL...		Metric Measurem...	Valid
Meta Data		urn:c...	Observed Value	90.0
VMD	MDC_DEV_ANALY_SAT...		Descriptor Version	0
Channel	MDC_DEV_ANALY_SAT...		Availability	CONTINUOUS
Numeric Metric	MDC_PULS_OXIM_SAT...	90.0	Metric Category	MEASUREMENT
Channel	MDC_DEV_PULS_CHAN		Resolution	1.0
Channel	MDC_DEV_PLETH_CHAN		Technical Ranges	0,00-127,00 (stepWidth...
Channel	MDCX_DEV_READINES...		Unit Code Id	MDC_DIM_PERCENT
SystemContext				

Our Tools

Device Modeler

- ▶ Model the device containment tree including **constraints for metrics**
- ▶ Generate network interface code including **runtime monitors**



Swiss Army Knife

- ▶ Debug devices in network and **manipulate them** through SCO
- ▶ Specify constraints for metrics and **execute monitors** checking those

Element	Term Code	Info	Attribute	Value
Network				
▶ NDS	MDC_DEV_ANALY_SAT...	UPFL...	Element	Numeric Metric
Raw MDDescription			Term Code	MDC_PULS_ORM_SAT...
Activate Operation D...	MDCX_START_IDENTF...		Handle	sp02
Activate Operation D...	MDCX_STOP_IDENTIF...		Observation Time	2018-05-02T23:00:55.301
Meta Data			Metric Measurement...	Valid
▶ VMD	MDC_DEV_ANALY_SAT...	UPFL...	Observed Value	95.0
▶ Channel	MDC_DEV_ANALY_SAT...		Descriptor Version	0
Numeric Metric	MDC_PULS_ORM_SAT...	MDL...	Availability	CONTINUOUS
▶ Channel	MDC_DEV_PULS_CHAN		Metric Category	MEASUREMENT
▶ Channel	MDC_DEV_PLETH_CHAN		Resolution	1.0
▶ Channel	MDCX_DEV_READINES...		Technical Ranges	0.00-127.00 [stepWidth...
▶ SystemContext			Unit Code Id	MDC_OIM_PERCENT

Specifying Temporal Properties

- ▶ Specify temporal behavior of metric's observed value, invocation state, ...
- ▶ Smart Assertion Logic for Temporal Logic (SALT, Bauer, Leucker, Streit 2006) ...
- ▶ Impartial Anticipation using LTL₃ Semantics (Bauer, Leucker, Schallhart 2011)

Example

```
assert "brightness.value < 70"  
until "readiness_state.value = 'READY'"
```

OSCP Device Modeler Workbench

OSCP Device Modeler XML Descriptions Medical Devices adm@tellinghusen.de

CREATE A NEW MEDICAL DEVICE

Workbench

Edit TestDevice

Zoom: 50%


The diagram illustrates a hierarchical tree structure for a medical device. The root node is 'TestDevice'. It branches into 'Device' and 'Component'. 'Device' further branches into 'Device_1' and 'Device_2'. 'Component' branches into 'Component_1' and 'Component_2'. Each node contains fields like Name, Number, Code, and Type. A vertical toolbar on the right contains buttons for VMD, CHN, MET, CP, CTX, ASY, ACO, ASG, and MCN. A home button is at the bottom right.

OSCP Device Modeler Workbench


OSCP Device Modeler XML_Descr

CREATE A NEW MEDICAL DEVICE Modeler Workbench


Edit TestDevice




OSCP Device Modeler




OSCP Device Architecture




OSCP Device Monitoring




OSCP Device Monitoring



OSCP Device Monitoring



OSCP Swiss Army Knife



UNIVERSITÄT ZU LÜBECK
INSTITUTE FOR SOFTWARE ENGINEERING
AND PROGRAMMING LANGUAGES

OSCP Swiss Army Knife

Manage Monitors


Show Monitor Events

Element	Term Code	Info	Attribute	Value
▼ Network			Element	Numeric Metric
▼ MDS	MDC_DEV_ANALY_SAT...	urn:u...	Term Code	MDC_PULS_OXIM_SAT...
Raw MDDescription			Handle	spo2
Activate Operation D...	MDCX_START_IDENTIF...		Observation Time	2016-05-02T23:00:55.301
Activate Operation D...	MDCX_STOP_IDENTIFI...		Metric Measurem...	Valid
Meta Data		urn:u...	Observed Value	90.0
▼ VMD	MDC_DEV_ANALY_SAT...		Descriptor Version	0
▼ Channel	MDC_DEV_ANALY_SAT...		Availability	CONTINUOUS
Numeric Metric	MDC_PULS_OXIM_SAT...	90.0	Metric Category	MEASUREMENT
▶ Channel	MDC_DEV_PULS_CHAN		Resolution	1.0
▶ Channel	MDC_DEV_PLETH_CHAN		Technical Ranges	0,00–127,00 (stepWidth...
▶ Channel	MDCX_DEV_READINES...		Unit Code Id	MDC_DIM_PERCENT
▶ SystemContext				

OSCP Swiss Army Knife

- ▶ Discover all active devices in network
- ▶ Show interface description
- ▶ Show current values
- ▶ Allow manipulation through service model
- ▶ **Allow attachment of monitors**

Monitor Injection in Swiss Army Knife



Demo

Demo

Conclusion

- ▶ IEEE 11073-SDC is upcoming standard for interconnection of medical point-of-care devices.
- ▶ Devices announce their interface to the local network.
- ▶ Monitors can be attached to devices' metrics.
- ▶ Monitors enforce own device and remote controlling device to satisfy the interface description.
- ▶ Swiss Army Knife discovers devices.
- ▶ Swiss Army Knife manipulates devices' states.